-  **DE Bedienungsanleitung**
-  **GB Instruction Manual**
-  **FR Mode d'emploi**
-  **IT Istruzioni per l'uso**
-  **ES Instrucciones de uso**

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General Information

About this Instruction Manual

Please read the safety instructions in this manual carefully. To avoid damage to the unit and the risk of injury, only use this product as described in the manual.

Keep the instruction manual handy so that you can easily look up information on all the functions.

DANGER!

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

CAREFUL!

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

NOTE!

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

Intended Use

This product is intended only for private use.

It was developed for the magnified display of things in nature.

General Warning

DANGER!

Tools with sharp edges and points are often used when working with this device. For this reason, store this device and all accessories and tools in a location that is out of the reach of children. There is a RISK OF INJURY!

DANGER!

This device contains electronic components which operate via a power source (power supply and/or batteries). Do not leave children unattended when using the device. Only use the device as described in the manual, otherwise you run the RISK of an ELECTRIC SHOCK.

DANGER!

Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short circuit, fires or even explosions.

DANGER!

Never bend, pinch, pull or run over the power and connecting cables or extensions/adapters. Protect the cables from sharp edges and heat. Before operation, check the device, cables and connections for damage.

Never use a damaged unit or a unit with damaged power cables, etc. Damaged parts must be exchanged by an authorized service centre immediately.

DANGER!

Children should only use the device under adult supervision. Keep packaging material (plastic bags, rubber bands, etc.) out of the reach of children! There is a RISK OF CHOKING!

Position your device so that it can be disconnected from the power supply at any time. The wall socket should always be located near the device and be easily accessible, since the plug on the power cord serves as a disconnecting device for the power supply.

CAREFUL!

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

NOTE!

Do not disassemble the device. In the event of a defect, please contact your dealer. He will contact the Service Centre and can send the device in to be repaired, if necessary.

Do not expose the device to temperatures above 40°C.

All parts (Fig. 1-7):

- 1 Eyepiece
- 2 Barlow lens 2x
- 3 Eyepiece holder
- 4 MikrOkular
- 5 Microscope head
- 6 Objective-revolver
- 7 Mikroskoptisch
- 8 Focus wheel
- 9 Filter wheel
- 10 Electric illumination
- 11 Brightness control
- 12 Microscope base
- 13 Slide Stage
- 14 Slide Stage retaining screw
- 15 Slide Stage adjustment (forwards/backwards)
- 16 Slide Stage adjustment (left/right)
- 17 Vernier scale
- 18 Stage clips
- 19 Specimen on a slide
- 20 On/off switch
- 21 LED status
- 22 Power supply cord

Accessories (Fig. 8-10):

- 1 MikrOkular
- 2 Software CD-ROM for MikrOkular
- 3 5x WF Eyepiece
- 4 10x WF Eyepiece
- 5 16x WF Eyepiece
- 6 Barlow lens 2x
- 7 Aluminium accessory box
- 8 Power supply
- 9 Specimen box with prepared slides, slides and cover glass

General/Location

Before you set up the microscope, you must choose a suitable location.

First you must make sure that your microscope is on a stable and solid surface.

An electricity supply is required for observation with the electric illuminator.

1. Stage Installation

The microscope's slide stage (Fig. 1, 13) is simple to attach. There are two small pegs on the bottom of the slide stage that fit directly into the holes on the microscope stage (Fig. 1, 7). The stage retaining screw (Fig. 1, 14) attaches the slide stage to the microscope stage.

2. Electronic Illumination

The Erudit is equipped with battery-operated LED lighting. This allows you to use your microscope without a power supply.

The microscope is operated with the mains adapter provided (Fig. 10, 8) to charge the battery or during use at home. To do so, connect the power cable to the microscope's power connection (Fig. 5, 22) and the adapter to a power point.

Press the on/off switch (Fig. 6, 20) to activate. The brightness of the lighting can be adjusted with the brightness control (Fig. 6, 11). The batteries begin charging after the microscope is connected to the power supply. The status LED turns red (Fig. 5, 21). When the battery is charged, the Status LED changes from red to green (Fig. 5, 21). Now you can use the microscope without the power supply.

The LED lighting lasts up to 48 hours when in battery operation. It can last longer, depending on the intensity of brightness used.

3. Microscope setup

The microscope's eyepiece (Fig. 1, 5) will now be prepared for the first observation. First, rotate the eyepiece into a convenient position.

Begin every observation with the lowest magnification. Place the microscope's table (Fig. 1, 7) with the focus wheel (Fig. 1, 8) into the lowest position and rotate the objective revolver (Fig. 1, 6) until it locks on the lowest magnification (4x)

Insert the 5x eyepiece (Fig. 9, 3) in the Barlow lens (Fig. 9, 6).

Take care, that the Barlow lens is inserted completely in the eyepiece holder (Fig. 1, 3).

4. Observation

After you have set up the microscope with the corresponding illumination, the following principles are important:

Begin each observation with the lowest magnification, so that centring and positioning as well as focussing of the object to be viewed is easier.

The higher the magnification, the more light is required for good picture quality.

Begin with a simple observation.

With the mechanical desk (Fig. 1, 13) you can position your specimen (Fig. 3, 19) smoothly and very precise.

The object is placed between the clamps on the mechanical desk.

Move the object, with help of the adjustment knobs, (Fig. 3, 15/16) to centre it under the objective.

With the built in vernier (Fig. 3, 17) at both axes you can now specifically set and shift the object. It can now be viewed with different magnifications.

Look through the eyepiece (Fig. 1, 1) and turn carefully the focusing wheel (Fig. 1, 8) until you can see a sharp picture.

Now you can get a higher magnification, while you slowly pull the Barlow lens (Fig. 4, 2) out of the eyepiece holder (Fig. 4, 3). If the Barlow lens is in the most outward position, the overall magnification of the microscope is multiplied by a factor of 1,6x.

For still higher magnification you can put the 16x eyepiece (Fig. 9, 5) into the eyepiece holder and set the objective revolver on higher magnification (10x / 60x). In this setting the overall magnification may also be multiplied by a factor of 1,6x by pulling out the Barlow lens.

Filter wheel:

The filter wheel (Fig. 1, 9) enables you to change the colour filter to improve the contrast of your specimens. The filter wheel has five different colour filters and one filter-free setting.

TIP:

Depending on the preparation higher magnifications do not always lead to better pictures.

Please notice:

When changing magnification (eyepiece or objective lens changes, pulling out of the Barlow lens) the sharpness of the image must be newly defined by turning the focusing wheel (Fig. 1, 8).

NOTE:

Please be very careful when doing this. When you move the mechanical plate upwards to fast the objective lens and the slide can touch and get damaged.

5. Viewed Object – condition and preparation

Condition

With an ordinary magnifying glass, we preferably view obscure items like little animals, plants etc. Here light shines on the viewed item, gets reflected and gets through the lens into our eye (Top illumination principle). With a biological microscope, a so called transmitted light microscope, only transparent items can be observed. Light gets through the item, becomes magnified by the objective and the eyepiece and gets in our eye. Many small water organisms, small plant parts (e. g. moss leaves) and finest animal components are naturally transparent; other ones must be accordingly prepared. By means of a pre-treatment or penetration with suitable materials (media) some specimens can be rendered transparent or we cut finest slices off of them (hand cut, Microtome cut) and then examine these. The following instructions make us familiar with these methods.

Preparing a thin cut of a specimen

As already before mentioned, very thin cuts have to be prepared from a non transparent object (plant rods, roots, etc.). In order to get satisfying results, we need some wax or paraffin to coat and harden the specimen if it is soft. If no such material is contained in your microscope set, then you can use a normal candle. The wax is given to a pot and melted over a flame. The object is then dipped several times into the liquid wax. Let the wax become hard. With a microtome or a knife / scalpel (caution!!!) thin slices are cut off of the object coated with wax. These cuts are put on a glass slide and covered with a cover glass.

Manufacture of an own preparation

Put the object which shall be observed on a glass slide and pipette a drop of distilled water on the object. Set a cover

glass (available in a well sorted hobby shop) perpendicularly at the edge of the water drop, so that the water runs along the cover glass edge. Lower now the cover glass slowly over the water drop.

DANGER!

Be extremely careful when using the knives/scalpels or the MicroCut. There is an increased risk of injury due to the sharp edges!

6. Experiments

If you made yourself familiar with the microscope already, you can accomplish the following experiments and observe the results under your microscope.

6.1 Newspaper print

Objects:

1. A small piece of paper from a newspaper with parts of a picture and some letters
 2. A similar piece of paper from an illustrated magazine
- Use your microscope at the lowest magnification and use the preparation of the daily paper. The letters seen are broken out, because the newspaper is printed on raw, low quality paper. Letters of the magazine appear smoother and more complete. The picture of the daily paper consists of many small points, which appear somewhat dirty. The pixels (raster points) of the magazine appear sharp.

6.2 Textile fibers

Items and accessories:

1. Threads of different textiles: Cotton, line, wool, silk, rayon, nylon etc..
2. Two needles

Each thread is put on a glass slide and frayed with the help of the two needles. The threads are wetted and covered with a cover glass. The microscope is adjusted to a low magnification. Cotton fibres are of vegetable origin and look like a flat, turned volume. The fibres are thicker and rounder at the edges than in the centre. Cotton fibres consist primary of long, collapsed tubes. Linen fibres are also of vegetable origin; they are round and run in straight lines. The fibres shine like silk and exhibit countless swellings at the fibre pipe. Silk is of animal origin and consists of solid fibres of smaller diameter contrary to the hollow vegetable fibres. Each fibre is smooth and even and has the appearance of a small glass rod. Wool fibres are also animal origin; the surface consists of structures which overlap like shingles, and appear broken and wavy. If it is possible, compare wool fibres of different weaving mills. Consider thereby the different appearance of the fibres. Experts can determine the country of origin of the wool by microscopic observation. Rayon (artificial silk) is manufactured by a long chemical process. All fibres show hard, dark lines on the smooth, shining surface. The fibres ripple themselves or crinkle after drying in the same condition. Determine common features and differences by examining many different fibres.

7. Setup MicrOcular

NOTE!

The MicrOcular only works without the Barlow lens supplied. Magnification setting is changed by using the MicrOcular and must be reset by focussing anew.

Remove the Barlow lens (Fig. 2, 2) with the currently used eyepiece from the tube (Fig. 2, 4) and put the MicrOcular (Fig. 4, 3) in place instead.

8. Installation and use of the software

NOTE!

Please do not yet connect MicrOcular and your PC. Please go through the items below in sequence.

8.1. Software and Installation Information

A software CD is included with your microscope. The software and drivers on this CD must be installed on your computer in order to use your MicrOcular. Once you connect the MicOcular to your computer after installation, you can view pictures on and save them to your computer using Photomizer SE. To install the software and driver correctly, simply follow the installation steps.

8.2. Installation der Software

1. Important: Before inserting the CD, first plug the USB cable into the USB port on your PC. Windows will now recognize the new device, and will indicate this in a notification window. Now please click on "Cancel", since we will have the driver for the device automatically installed later.
2. Now insert the CD-ROM that came in your package into the CD/DVD drive of your computer. The installation menu starts automatically. If it does not, go to the Windows explorer and select the CD/DVD drive (most of the time, it's the "D" drive, but it can have another letter). From there, start the file "menu.exe" by double-clicking with the left mouse button.

8.2.1. Driver installation

To install the driver software, click on the menu point "Install Driver" with the mouse cursor. Then follow the installation program instructions.

During the software installation the correct driver for your operating system will be installed automatically. No manual input is needed. In rare cases the device may not be recognised by your computer. As a rule you need then only install the driver manually from the CD. If this fails please refer to the troubleshooting chapter that follows.

8.2.2. Installation of the Image Editing Software Photomizer SE

The image editing software "Photomizer SE" is located on the software CD. You can edit your pictures here.

1. To install the software, click once on the menu point "Install Software" with the left mouse button.
 2. The Photomizer Software requires Microsoft .NET Framework 3.5, which it will install if it is not already on your system. If it is already on your system, you may skip to step 5.
 3. In the welcome window, please accept the Microsoft license agreement, and then click "Install". The installation can take a few minutes.
 4. As soon as everything is installed, click "Finish".
 5. Now you will be presented with a choice, in which you can choose your language.
- Make your selection and confirm it by clicking „OK“.
6. When you see "Welcome", click on „Next“.
 7. In the next window, you will be asked for the „Destination Folder“. Here, just click on „Next“.
 8. Now, the window with the Setup-Status will appear – here, a progress bar will inform you about the ongoing installation. This process can take a few minutes.
 9. The window "Photomizer is being installed" appears. Click on "Finish". The installation ends.

TIP:

To use the MicrOcular long-term we recommend it always be connected to the same USB port.

8.3. Troubleshooting: the device isn't recognised by the computer (USB)

1. This occurs in rare cases with some OEM Windows Vista or XP versions. The device is not defective. Updating the notebook or PC USB driver usually suffices.
2. If not, proceed as below. (The maker cannot accept any liability whatsoever in these cases¹).

Delete ghost devices. These are devices not actually currently connected to your computer. The reason for this is that Windows generates a system settings entry for each USB device (e.g. USB Stick) on each port. It then seeks that/those device/s each time it's started. This slows down USB device recognition and Windows booting and causes non-recognition of USB ports.

Cleaning up the system environment and deleting devices that are no longer in use helps here. To do so right click WORKPLACE (COMPUTER in Vista) and open the PROPERTIES menu. Then click the ADVANCED tab (ADVANCED SYSTEM SETTINGS in Vista) and then ENVIRONMENT VARIABLES. Choose the NEW option in ENVIRONMENT VARIABLES. Enter the following under VARIABLE NAME.

devmgr_show_nonpresent_devices

Under VARIABLE VALUE enter „1“.

Confirm your input with OK and then reboot the computer.

Change to the device manager. In the VIEW menu activate the SHOW DEVICES OPTION. The ghost devices previously not shown will now be displayed in light grey. Check the various categories including USB, memory size, etc. Only delete entries for devices you no longer use.¹

¹Source: Microsoft Knowledge Base

<http://support.microsoft.com/kb/315539/en>

9. Using the MicrOkular

9.1. Preparation

1. Slide a specimen under your microscope and focus on it.
2. Remove the eyepiece and Barlow lens from the eyepiece support and the dust cap from your MikrOkular and install same in the eyepiece supports instead of the Barlow lens.
3. Start your PC if you haven't yet and connect your MikrOkular to the USB port of your computer.

9.2. Showing and saving MicrOkular images on your PC

1. Start the Photomizer SE Software.
2. Click „Open camera“
3. In the event that you have connected more than one unit, you can choose the desired unit in the subsequent selection. Click on „USB 2.0 Webcam“. If only one device is connected, this step is omitted.
4. The camera image should now be visible on your screen. Focus the microscope image.
5. Click „Capture“ to record an image. It will then be shown on the right.
6. Click on it to select it and then click „Transfer image“.
7. This takes you to the Photomizer SE software.
8. File - Save as

9.3. The Photomizer SE Software

If you need help in the „Photomizer SE“ programme click „?“ and then „Open help“. If you need further assistance please visit the maker's home page at www.photomizer.net

NOTES on cleaning

Separate the device from the power supply before cleaning (remove plug).

Only use a dry cloth to clean the exterior of the device.

NOTE!

Do not use any cleaning fluid to avoid damaging the electronics.

Clean the lenses (eyepiece and/or lens) with a soft and lint-free cloth only (e.g. microfibre).

NOTE!

Do not apply excess pressure to the cloth so as to avoid scratching the lenses.

To remove more stubborn dirt, moisten the cleaning cloth with an eyeglass-cleaning solution and wipe the lenses gently.

Protect the device from dust and moisture! After use - in particular in situations of high humidity - let the device acclimatize for a short period of time, so that the residual moisture can dissipate. Remove the dust cover and store it in the included bag.

Troubleshooting

Error

No picture visible

Solution

- switch light on
- put condenser lens in place
- adjust focus

Picture flickers while viewing with PC eyepiece

- if necessary, adjust resolution of the video graphics board

Software installation reports "not XP approved"

- confirm by clicking "OK"

Technical Data

System requirements for PC eyepiece

Minimum system requirements: PC with an Intel Pentium IV processor or higher; Windows XP with Service Pack 3 (on CD-ROM), Windows Vista or Windows 7 - with DirectX 9.x (on CD-ROM), .NET Framework 3.5 (on CD-ROM); min. 1024 MB RAM; min. 500 MB free hard drive space; free USB port.

Microscope

Monocular head

Mechanical desk (with vernier scale)

Eyepieces 5x, 10x, 16x

Illumination: Accu, LED

Magnification table

Eyeiece	Objective	Magnification	(with Barlow lens)
5x	4x	20x	40x
5x	10x	50x	100x
5x	40x	200x	400x
10x	4x	40x	80x
10x	10x	100x	200x
10x	40x	400x	800x
16x	4x	64x	128x
16x	10x	160x	320x
16x	40x	640x	1280x

Technical alterations and errors excepted.

Disposal

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage! As per the Directive 2002/96/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner. Empty old batteries must be disposed of at battery collection points by the consumer. You can find out more information about the disposal of devices or batteries produced after 01.06.2006 from your local waste disposal service or environmental authority.

Please take the current legal regulations into account when disposing of your device. You can get more information on the proper disposal from your local waste disposal service or environmental authority.

Declaration of Conformity

Product Type: Transmission-type microscope
Product Name: BRESSER Erudit MO
Article No.: 51-10000

Bresser GmbH has formed a „Conformity Declaration“ in line with the applicable guidelines. This declaration is available for review.

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Warranty

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Microscopes which become defective during the warranty period can be returned to the dealer where the microscope was bought. The repaired microscopes or a new pair will then be returned to you, free of charge. In the case of defects which occur after the end of the warranty period, the microscopes can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the microscopes carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This guarantee does not imply any restriction of your statutory rights.

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