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BRESSER Science ADL-601P
Reflection-type and
Transmission-type Microscope
with Polarization Attachment

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CAUTION!

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

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Do not dispose of electric equipment together with household waste material! In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of its life must be collected separately and returned to an environmentally compatible recycling facility. Discharged batteries and damaged re-chargable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority.

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General Information

About this Instruction Manual

Please read the safety instructions in this manual carefully. To avoid damage to the unit and the risk of injury, only use this product as described in the manual.

Keep the instruction manual handy so that you can easily look up information on all the functions.

**DANGER!**

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

**CAREFUL!**

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

**NOTE!**

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

**CAREFUL!**

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

**NOTE!**

Do not disassemble the device. In the event of a defect, please contact the Service Centre in your country (see the warranty card).

Do not expose the device to temperatures above 45°C.

Intended Use

This product is intended only for private use.

It was developed for the magnified display of things in nature.

**General Warning****DANGER!**

Tools with sharp edges and points are often used when working with this device. For this reason, store this device and all accessories and tools in a location that is out of the reach of children. There is a RISK OF INJURY!

**DANGER!**

This device contains electronic components which operate via a power source (power supply and/or batteries). Do not leave children unattended when using the device. Only use the device as described in the manual, otherwise you run the RISK of an ELECTRIC SHOCK.

**DANGER!**

Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short circuit, fires or even explosions.

**DANGER!**

Never bend, pinch, pull or run over the power and connecting cables or extensions/adapters. Protect the cables from sharp edges and heat. Before operation, check the device, cables and connections for damage.

Never use a damaged unit or a unit with damaged power cables, etc. Damaged parts must be exchanged by an authorized service centre immediately.

**DANGER!**

Children should only use the device under adult supervision. Keep packaging material (plastic bags, rubber bands, etc.) out of the reach of children! There is a RISK OF CHOKING!

The “BRESSER Science ADL-601P” is equipped with reflected and transmitted light systems and a polarization attachment. With the planachromatic objectives and wide field eyepieces, the observer can get clear images of a wide field

of view. The instrument is suitable for knowing something about the internal construction of materials and for condition monitoring of different mechanical apparatus in the colleges and the factory.

Table of Contents

- I. Components
- II. Technical Data
- III. Installation
- IV. Observing operation
- V. Maintenance, Declaration of conformity, Warranty & Service

I. COMPONENTS

Fig. 1

- 1 Binocular tube
- 2 Photo tube
- 3 Polarizer (incident illumination)
- 4 Color filter wheel
- 5 Holding screw for lamp housing
- 6 Condensor adjusting lever (incident illumination)
- 7 Lamp housing (incident illumination)
- 8 Adjusting knob incident light bulb (tilt)
- 9 Adjusting knob incident light bulb (up/down)
- 10 Cross table adjustment (cross moving)
- 11 Cross table adjustment (lengthwise)
- 12 Reflected illumination dimmer
- 13 Base
- 14 Head holding screw
- 15 Analyzer adjusting rod

Fig. 2

- 16 Limit knob (up stop)
- 17 Power supply plug
- 18 AC adapter
- 19 Socket for incident light illumination
- 20 Knob for tensional adjustment of focusing
- 21 Coarse focus
- 22 Fine focus

Fig. 3

- 23 Aperture diaphragm adjustment lever (incident illumination)
- 24 Field diaphragm adjustment lever (incident illumination)
- 25 Switch for photo mode / observation
- 26 Analyzer rod
- 27 Nosepiece
- 28 Objective
- 29 Slide holder
- 30 Stage
- 31 Transmitted illumination dimmer
- 32 On/Off switch
- 33 Condensor up / down knob
- 34 Aperture diaphragm adjustment lever (transmitted illumination)

Fig. 4

- 35 Condenser centering screw
- 36 Collector
- 37 Polarizer (transmitted illumination)
- 38 Condenser holding screw
- 39 Filter holder
- 40 Aperture diaphragm adjustment lever (transmitted light)

II. TECHNICAL DATA

Objective

Eyepiece type	Magnification	Field of view diameter / mm	Focal length f / mm	Remarks
Wide field eyepiece (WF)	10x	18	25	
WF eyepiece with crosshair	10x			optional

Objective

Objective type	Magnification	Numerical aperture NA	Working distance WD / mm	Remarks
Planachromatic objective (no cover glass, except for 100x)	4x	0.10	17.9	optional
	5x	0.12	18.3	
	10x	0.25	8.9	
	40x	0.60		
	60x	0.75		
	100x	1.25 (oil)		optional

Total Magnification

Objectives:	4x	5x	10x	40x	60x	100x
Total Magnification						
Eyepiece:						
10x	40x	50x	100x	400x	600x	1000x
16x	64x	80x	160x	640x	960x	1600x

- Illumination: 12 V, 20 W halogen lamp, brightness adjustable.
- Power supply: AC 230 V / 50 Hz.
- The included filters of the reflected light system are: polarized, frosted glass, blue, green and red.

- Stage movement range: 50 mm x 75 mm.
- Minimum division of fine focusing scale: 0.002 mm.
- Adjustment range of interpupillary distance: 53-75 mm.
- Anti-fungus: Binocular tube and microscope head.


Fig. 5

III. INSTALLATION

1. Incident illumination (Fig. 5)

Take off the protective cap of the auxiliary lens under the reflected light illumination tube, install the reflected light illumination tube on the socket (3) of the body (5), and tighten up the tube holding screw (4).

2. Lamphouse (Fig. 5)

Loosen the lamphouse holding screw (1) and install the lamphouse (2) into the end of the reflected light tube. Then tighten the lamphouse holding screw (1) again.

IV. OBSERVING OPERATION

1. Incident illumination (Fig. 1, 3 & 6)

1. Plug the power supply plug into the microscope (19).
2. Push the power switch (32) to "I"; then the light power is turned on.
3. Unscrew the 10x objective (28) and turn the nosepiece (27) to its anchor point, so that the empty hole of the 10x objective is luminous and in the working position.
4. Put a piece of white paper on the stage (30). You should see an image of the luminous hole on the paper.
5. Adjust the lamp with the respective adjusting knobs (8, 9), until the image of the light bulb filament is centered on the bright light spot on the paper.
6. If the image of the filament is not clear, you should adjust the condenser adjusting lever (6) to make the image clear.
7. Install the 10x objective and put an object on the stage.
8. Adjust the coarse/fine focusing knob (21 & 22) in order to make the image of the object clear.
9. If the brightness of the field is not equal, you can move the condenser adjusting lever (6) slightly to make the brightness of the field equal.
10. Turn the field diaphragm (24) to minimum; adjust the field diaphragm centering screw (34) to superposition the center of the field diaphragm and the center of the field. Then open the field diaphragm until you cannot see the image of the field diaphragm in the field.
11. When you use other objectives with another magnification, you can adjust the aperture diaphragm and the brightness of illumination to obtain a satisfactory contrast.
12. In order to obtain an image of good quality, you can use a piece of frosted glass or filters (4) according to your need.


Fig. 6

2. Transmitted illumination (Fig. 4 & 7)

Fig. 7


- 1 Switch for observation / photo mode
- 2 Analyzer adjusting lever

1. Push the power switch (32) to "I"; then the illumination power is turned on.
2. Adjust the transmitted illumination dimmer (31) and the aperture diaphragm adjustment lever (40), turn the condenser up/down knob (33) to move the condenser.
3. If applicable, move the observation / photo switch (25) into the observing position and push the analyzer with the analyzer rod (26) into the light path.

3. Interpupillary distance adjustment (Fig. 8)

Fig. 8

- 1 Binocular Tube
- 2 Diopter adjustment ring


Put a specimen on the stage and focus to obtain the image of it. Adjust the interpupillary distance of the binocular tube (2) by "folding" until the right and left field of view can be composed to one.

4. Diopter adjustment (Fig. 8)

Turn the 40x objective to the observation position. Firstly, observe through the right eyepiece with the right eye; adjust the coarse and fine focusing knobs to image clearly. Secondly, observe through the left eyepiece with the left eye; adjust the diopter adjustment ring (1) to image clearly.

5. Coarse / Fine focusing (Fig. 9)

Fig. 9

- 1 Fine focusing knob
- 2 Coarse focusing knob
- 3 Knob for tensional adjustment of focusing
- 4 Limit knob (up stop)

The instrument uses a coarse/fine focusing mechanism. The knob for tensional adjustment of focusing (3) is used for adjusting the tension of the coarse focusing knob (2) to prevent the stage from naturally sliding down and to adjust the knob to be comfortable. The limit knob (4) prevents accidental specimen/objective contact. After locking the limit knob (by turning up and fastening) at a chosen height of the stage, you cannot focus the object closer to the objective using the coarse focusing knob. So, the object is protected against damage. (1) is fine focusing knob.

6. Incident illumination observation system (Fig. 1, 3 & 4)

Put an object on the stage (30), select the 10x objective and observe with it, turn the coarse and fine focusing knobs (21 & 22) to make the image of the object clear. Turn the field diaphragm down by turning the field diaphragm adjustment lever (24), then adjust the field diaphragm centering screw (34) to make the field diaphragm center coinciding with the optical axis. Then you can turn the field diaphragm adjustment lever (24) to make the field diaphragm bigger than the field of view. Turn the filters turn-plate (4) to select a filter.

7. Transmitted illumination observation system (Fig. 3 & 4)

Put a preparation on the stage (30), select the 10x objective and observe with it, turn the coarse and fine focusing knobs (21 & 22) to make the image of the object clear. When you observe with an objective of different magnification, you should focus by the fine focusing knob (22). Then put the filter into the filter holder (39).

8. Observation with polarized illumination (Fig. 1, 3 & 4)

1. Push the analyzer rod (26) to bring the analyzer into the working position.
2. When you observe with reflected illumination, put the polarizer (3) into its insertion slot and turn the polarizer (3). You can observe the variant light contrast in the field. When the polarized light is in quadrature (polarizer and analyzer are orthogonal), you can observe the specimen on the stage with polarized light.
3. When you observe with transmitted illumination, you should put the polarizer onto the collector (36) and turn the polarizer to observe the variant light contrast in the field until the polarized light is in quadrature.

9. Exchange of the lamp

9a. Incident illumination (Fig. 10)


Fig. 10

- 1 Lamp
- 2 Lamp holder
- 3 Lamp holder fixing screw

1. Switch off the power supply and put out the plug from mains socket.
2. Loosen the lamp holder fixing screw (3) and open the lamp housing. Put out the lamp holder (1). Be careful, as the lamp may be hot!
3. Pull out the old lamp (1) from the lamp holder (2) and install a new lamp. When you install the new lamp, be careful not to touch the glass with your fingers. The new lamp should be supplied in a plastic protective envelope. If not, use a tissue or another medium to grasp the lamp. This will prevent contamination from your hand, which can reduce the lamp's intensity and life.
4. Install the lamp holder (2), plug in the plug and switch on the power supply again. Then adjust according to paragraph 1 ("Incident illumination").

9b. Transmitted illumination (Fig. 11)


Fig. 11

- 1 Microscope base (bottom)
- 2 Lamp socket
- 3 Lamp
- 4 Lamp centering screw
- 5 Lamp door retaining screw

1. Switch off the power supply and put out the plug from mains socket.
2. Incline the microscope, loosen the lamp door retaining screw (5) on the middle part of the bottom (1) and open the lamp door.
3. Pull out the old lamp (3) from the lamp socket (2). Be careful, as the lamp may be hot!
4. Insert the new lamp into the lamp socket (2). Notice the properly touching; take care not to touch the lamp with bare fingers. E. g., use the protective envelope of the lamp or a tissue, in order to grasp the bulb.
5. Reinstall the lamp door with lamp stand board on the bottom by tightening the screw (6).
6. Plug in the plug, switch on the power supply again and observe. If the illumination is not even, you can loosen the lamp centering screw (4) and adjust the position of the lamp.

V. MAINTENANCE

1. Sweep the lens
Sweep the lens by lens tissue or soft fabric dampened with alcohol.
2. Clean the painted parts
The dust on the painted parts can be removed by gauze. For the grease spots, the gauze immersed slightly with alcohol is recommended. Do not use organic solvents such as ether or other thinner etc. for cleaning the painted parts or plastic components.
3. Avoid disassembling the microscope
Because of being a precise instrument, do not disassemble the microscope casually. That may cause serious damage to its performance.
4. When not in use
Cover the microscope with the dust cover and place it there, where it is dry and mouldless. We suggest the storage of all objectives and eyepieces in a closed container with drying agent.

EEC CONFORMITY EXPLANATION

Bresser GmbH, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 61326: 1997
EN 61000-3-2

Product description: Reflection-type and Transmission-type
Microscope with Polarization Attachment
Model: BRESSER Science ADL-601P

Rhede, 20. Juli 2007

Bresser GmbH

Helmut Ebbert
Managing director

SERVICE & WARRANTY

The regular guarantee period is 2 years and begins on the day of purchase. To benefit from an extended voluntary guarantee period as stated on the gift box, registration on our website is required.

You can consult the full guarantee terms as well as information on extending the guarantee period and details of our services at www.bresser.de/warranty_terms.

Would you like detailed instructions for this product in a particular language? Then visit our website via the link below (QR code) for available versions.

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Alternatively you can also send an email to manuals@bresser.de or leave a message on +49 (0) 28 72 – 80 74-220*. Please always state your name, precise address, a valid phone number and email address, as well as the article number and name.

*Number charged at local rates in Germany (the amount you will be charged per phone call will depend on the tariff of your phone provider); calls from abroad will involve higher costs.

DISPOSAL

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage! As per the Directive 2002/96/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner. Empty old batteries must be disposed of at battery collection points by the consumer. You can find out more information about the disposal of devices or batteries produced after 01.06.2006 from your local waste disposal service or environmental authority.

Please take the current legal regulations into account when disposing of your device. You can get more information on the proper disposal from your local waste disposal service or environmental authority.

Bresser GmbH

Gutenbergstr. 2 · DE-46414 Rhede · Germany

Tel. +49 (0) 28 72 - 80 74-210

Fax +49 (0) 28 72 - 80 74-222

www.bresser.de · service@bresser.de